

Jun 11, 2026

WATERS SELECT SERIES MRT MALDI Spot Analysis for Zooarchaeology by Mass Spectrometry (ZooMS) of collagenous material


DOI


<https://dx.doi.org/10.17504/protocols.io.n2bvjkwbnk5/v1>

Frankie Tait¹, Nicholas Michael², Karen Ruebens¹

¹Archaeological Proteomics Laboratory, Department of Archaeology, University of Reading, UK;

²Chemical Analysis Facility, University of Reading, UK

 **Frankie Tait**
University of Reading



OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvjkwbnk5/v1>

Protocol Citation: Frankie Tait, Nicholas Michael, Karen Ruebens 2026. WATERS SELECT SERIES MRT MALDI Spot Analysis for Zooarchaeology by Mass Spectrometry (ZooMS) of collagenous material. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n2bvjkwbnk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 21, 2026

Last Modified: June 11, 2026

Protocol Integer ID: 315439

Keywords: ZooMS, Palaeoproteomics, Mass Spectrometry, MALDI-TOF, collagen fingerprinting, Q-ToF, zooarchaeology by mass spectrometry, peptide mass fingerprinting of collagen type, peptide mass fingerprinting, select series mrt maldi spot analysis, peptide peak, mass spectrometry, waters mrt, maldi spot analysis capability, zooarchaeology, taxonomic identification, collagenous material this protocol, mrt, collagen type

Funders Acknowledgements:

UKRI
Grant ID: EP/Y037448/1

Abstract

This protocol gives guidance on how to set up and run a Waters™ Select Series Multi-Reflecting Time-of-Flight (MRT) with a MALDI source attachment for the purposes of Zooarchaeology by Mass Spectrometry (ZooMS) of collagenous material. ZooMS allows for taxonomic identifications to be made based on peptide mass fingerprinting of collagen type 1. The Waters MRT is an ultra-high resolution imaging instrument with MALDI spot analysis capabilities, allowing for high precision ZooMS to be carried out. This instrument has the potential to achieve sub-ppm mass accuracy, therefore, results can give peptide peaks accurate to at least one decimal place allowing for contaminants and non-diagnostic peptide markers to be distinguished.

Equipment	
Select Series MRT	NAME
MALDI	TYPE
Waters	BRAND
https://www.waters.com/nextgen/us/en/products/mass-spectrometry/mass-spectrometry-systems/select-series-mrt.html?srsId=AfmBOopwbKjMBGQ3IWwtUfEUhpYa8cC7Pwq8IRwSnPpKqvW61xH_cnuz#:~:text=Maximize%20MS%20imaging%20capabilities%20with%20DESI%20and%20MALDI&te	LINK
355 nm Nd:YAG laser, tightly focused laser beam profile; Quadrupole Multi-Reflecting Time of Flight (Q-ToF)	SPECIFICATIONS

Attachments



FILE

[_extern.inf](#)

12KB



PNG

[MALDI_Values.png](#)

259KB



PNG

[StepWave_Values.png](#)

134KB



PNG

[MSProfile_Values.png](#)

189KB



PNG

[IonGuides_Values.png](#)

197KB



PNG

[Transfer_Values.png](#)

155KB



PNG

[RF_Values.png](#)

107KB



FILE

[Sample_A.ms](#)

2.3MB



FILE

[Sample_B.ms](#)




2.2MB




Guidelines





1. See "Materials" for reagent recipes
2. Mass Spectrometry is sensitive to contaminants eg. detergents, plastics etc.
3. This protocol can be performed in a standard wet chemistry laboratory.
4. Wear relevant PPE as designated by your local laboratory.
5. Follow local lab guidelines regarding sample handling and storage.
6. Follow country and facility specific guidelines regarding the disposal of chemical waste.
7. This protocol is intended to give guidance to mass spectrometry and/or Waters technicians - instrument settings should not be altered without specialist training.

Materials





Reagents

0.1% TFA/ 5% Acetonitrile (v/v) in HPLC-MS water: Add  940 μL HPLC-MS water to a 1.5 mL microcentrifuge tube. Add  50 μL ACN. Carefully add  10 μL 10% TFA. Secure cap and vortex briefly.

Conditioning Solution; 0.1% TFA in 50:50 Acetonitrile and HPLC-MS water: for 500ml, add  250 mL ACN to  250 mL LCMS water and add  500 μL 100% TFA.

 α -Cyano-4-hydroxycinnamic acid (CHCA) | **VWR Catalog #SIALC2020-10G** :  10 mg/mL in Conditioning Solution. Ensure powder is fully dissolved either through vortexing or heating (or sonicating if using a glass vessel) - hold up to the light to be sure. Store at  Room temperature for no more than  24:00:00 . Store in the dark.

 Red Phosphorous (RedP) | **Merck MilliporeSigma (Sigma-Aldrich) Catalog #343242-5G**  25 mg/mL finely ground suspension in ACN.

 Waters™ MassPREP ADH Digestion Standard | **Waters Catalog #186002328** Resuspend vial in  1 mL of 0.1% TFA/ 5% ACN (v/v) in LCMS water. This will give a concentration of  1 pMol/ μL . Aliquots should be stored  -20 °C and only thawed immediately before use, do not continuously freeze-thaw.

 Methanol (MeOH) | **Fisher Scientific UK Catalog #10653963**

 Water | **Fisher Scientific UK Catalog #10777404**

 Acetonitrile (ACN) | **Fisher Scientific UK Catalog #10616653**

Equipment	
MALDI 480 Target Plate	NAME
Waters	BRAND
405018182	SKU
Polished Stainless Steel	SPECIFICATIONS

Equipment	
Select Series MRT	NAME
MALDI	TYPE
Waters	BRAND
https://www.waters.com/nextgen/us/en/products/mass-spectrometry/mass-spectrometry-systems/select-series-mrt.html?srltid=AfmBOopwbKjMBGQ3IWwtUfEUhpYa8cC7Pwq8IRwSnPpKqvW61xH_cnuz#:~:text=Maximize%20MS%20imaging%20capabilities%20with%20DESI%20and%20MALDI&te	LINK
355 nm Nd:YAG laser, tightly focused laser beam profile; Quadrupole Multi-Reflecting Time of Flight (Q-ToF)	SPECIFICATIONS


Protocol materials

 Water | Fisher Scientific UK Catalog #10777404

 Methanol (MeOH) | Fisher Scientific UK Catalog #10653963


 Red Phosphorous (RedP) | Merck MilliporeSigma (Sigma-Aldrich) Catalog #343242-5G

 Waters™ MassPREP ADH Digestion Standard | Waters Catalog #186002328

 α -Cyano-4-hydroxycinnamic acid (CHCA) | VWR Catalog #SIALC2020-10G

 Acetonitrile (ACN) | Fisher Scientific UK Catalog #10616653

Safety warnings

-  *It is recommended that for MALDI analysis on the MRT, 50 volts is added to the detector voltage that has been obtained via the automatic detector gain setup. This should be done in consultation with the manufacturer.


Before start

1. We recommend using dedicated glassware where possible for preparation and storage of reagents.
2. Peptide extracts for troubleshooting this protocol were produced using an Acid Insoluble ZooMS protocol modified from Buckley et al. 2009 and Welker et al. 2015.
3. We do not recommend spotting the peptide extracts more than 24 hours in advance.
4. Ensure the MALDI laser window, mirror and octopole are clean.


Spotting


1 Spot  Sample 1:1 with 

 α -Cyano-4-hydroxycinnamic acid (CHCA) | **VWR Catalog #SIALC2020-10G**


matrix (see Materials). Spot peptides first, then mix with CHCA by pipetting in a circular motion on the plate  1 μ L of each


Note

Final volume in each target well should be  2 μ L to cover full surface area

2 Spot the  1 pMol/ μ L aliquots of

 Waters™ MassPREP ADH Digestion Standard | **Waters Catalog #186002328** (see


Materials) 1:1 with CHCA in the same manner as the peptide extracts. Spot 8 replicates. Each spot has a final concentration of  500 fMol/ μ L

3 Allow the plate to dry on a heater block  35-50 °C to encourage crystallisation. Remove as soon as the last spot has dried.



4 The plate should be kept at room temperature in the dark if not running immediately.

Instrument Calibration

5 Ensure 

 Red Phosphorous (RedP) | **Merck MilliporeSigma (Sigma-Aldrich) Catalog #343242-5G**

(see Materials) is resuspended before and during this step by shaking gently.

6 Spot  1 μ L of RedP onto a calibrant spot 4-5 times, allowing  00:00:20 in between each dispense for the layer to dry.

7 Load the plate into the instrument and ensure the vacuum has been established.



8 Acquire data for the calibration using the following settings

Note

- Positive ; MRT ; Standard ; 1.0s ; 9600
- API Gas flow rate 350 mL/min
- Laser Rate 2500 Hz
- Laser Intensity 300
- Lens Focus 4.00 mm
- Step Size: X 1000 μm and Y 100 μm
- Scan Time 1.0 sec
- Decoding Threshold Automatic
- Function MS

9 Only accept the calibration if the Mean Prediction Error is ≤ 0.2 ppm.



Instrument Standard and Laser Optimisation

10 Test three laser intensities: 250, 300 and 350 to confirm optimum settings using the spotted ADH.

Note

The laser intensity needs to be optimised depending on variable environmental factors and laser wear and tear. This will also act as an instrument standard to ensure functionality before and after processing samples.

11 Acquire data using the following recommended settings



Note

- Positive ; MRT ; Standard ; API Gas ; 1.0s ; 9600
- Laser Rate 2500 Hz
- Lens Focus 4.00 mm
- Well Diameter: 2.15 mm ; Step Size: X 500 μ m and Y 150 μ m
- Scan Time 1.0 sec
- m/z range 1000 - 3500
- Decoding Threshold Low
- Function MS

Manual Quad Profile:


- Dwell time 20% at 1000 m/z
- Dwell time 40% at 1500 m/z
- Top end mass 3000 m/z
- Ramp time 20%

The RF controls below were applied to allow for higher peptide markers to pass through the quadrupole:

MALDI	StepWave	Quad/MS Profile/DRE	Ion Guides	Transfer	RF
RF Settings					
StepWave RF (V):	250	250			
Ion Guide RF (V):	700	700			
Trap RF (V):	200	200			
Driftcell 1 RF (V):	200	200			
Driftcell 2 RF (V):	200	200			
Transfer RF (V):	800	800			
Transfer RF Gain:	5				
RF Ramp Settings					
Ion Guide RF Ramp Enable:	OFF				
MALDI RF					
Hex RF (V):	350	349			
Hex RF Gain (V):	0				

See attachments for full settings and an example run file.

12




ADH Spectrum

Expected result

Key peptide peaks that should be observed are 2744 m/z and 3369 m/z to ensure higher mass peptides are being appropriately measured.

Sample Data Collection

- 13 Use the same settings as for the ADH instrument standard. Set the laser intensity to the optimised setting based on spectral quality of signal-to-noise and TOF MS ES+ counts.



Sample Spectrum

- 14 Ensure you set the instrument to acquire all sample data, followed by one ADH spot to confirm consistency throughout the run.
- 15 Each spot receives 30 scans which will be summed.

Data processing

- 16 The Waters raw data files (.raw) can be converted to mzML format using MSConvert from ProteoWizard (Chambers et al., 2012)

- 17 The individual scans can be summed into .txt files using an in-house developed Python script (MZSummer), which we make publicly available here:

<https://github.com/nm303909/MZSummer/tree/main>

The screenshot shows the MZSummer web interface with the following settings:

- Input Directory:** [Empty text box] **Browse**
- Output Directory:** [Empty text box] **Browse**
- **Sub-Scan Selector****
 - Base Peak Intensity Threshold: 5000
- **Sub-Scan Filter****
 - m/z Range to Process (Min, Max): 1000 [Empty] 3500
 - Filter by: Sub-Scan Filter Intensity Sub-Scan Filter S/N Threshold
 - Sub-Scan Filter Intensity: 5000
 - Sub-Scan Filter S/N Threshold: 3
- **Output Options****
 - Minimum Intensity of m/z Output: 100
 - Decimal Places for m/z values: 2
 - Include Headers in Output
- Buttons:** About, Help, Submit
- **Lock Mass Correction****
 - Enable Lock Mass Correction
 - Lock Mass m/z: 609.28066
 - Lock Mass Window (ppm): 20.0
 - Lock Mass Minimum Intensity: 10000

MZSummer user interface

See our standard settings in the screenshot above - these settings can be adjusted to suit your needs. (We recommend running this programme overnight as it can take several hours if running a full plate on a laptop)




- 18 After summing, any technical replicates can be merged in R using the MALDIquant and MALDIquantForeign packages and published scripts (Gibb and Strimmer, 2012; Mylopotamitaki et al., 2024)
- 19 Once you have fully processed the raw data into a workable file, we advise following this separate guide for completing spectral identifications:



<https://www.protocols.io/view/identifying-zooms-spectra-mammals-using-mmass-kqdg36rppg25/v1>. Some modifications may be necessary due to the ultra-high resolution of the spectra (e.g. you may consider unselecting the deisotoping panel for peak picking due to the narrow peaks).

See attachments for example summed and merged MRT spectra: Sample_A and Sample_B

Target Plate Cleaning

- 20 Rinse the plate  Methanol (MeOH) | Fisher Scientific UK Catalog #10653963 .
Visibly check the plate and use a cotton bud, damp with MeOH, to manually wipe any residue that was not easily rinsed away.
- 21 Sonicate the plate in 50:50 MeOH and  Water | Fisher Scientific UK Catalog #10777404 for  00:15:00
- 22 Dry in the fume hood either under regular extraction or with a Nitrogen gas tube.

Protocol references

- Buckley, M., Collins, M., Thomas-Oates, J., & Wilson, J. C. (2009). Species identification by analysis of bone collagen using matrix-assisted laser desorption/ionisation time-of-flight mass spectrometry. *Rapid Communications in Mass Spectrometry: RCM*, 23(23), 3843–3854. <https://doi.org/10.1002/rcm.4316>
- Chambers, M. C., Maclean, B., Burke, R., Amodei, D., Ruderman, D. L., Neumann, S., Gatto, L., Fischer, B., Pratt, B., Egertson, J., Hoff, K., Kessner, D., Tasman, N., Shulman, N., Frewen, B., Baker, T. A., Brusniak, M.-Y., Paulse, C., Creasy, D., ... Mallick, P. (2012). A cross-platform toolkit for mass spectrometry and proteomics. *Nature Biotechnology*, 30(10), 918–920.
- Gibb, S. and Strimmer, K., 2012. MALDIquant: a versatile R package for the analysis of mass spectrometry data. *Bioinformatics*, 28(17), pp.2270-2271.
- Mylopotamitaki, D., Fewlass, H., Zavala, E.I., Rougier, H., Sümer, A., Hajdinjak, M., Smith, G.M., Ruebens, K., Sinet-Mathiot, V., Xia, H., Hansen, J., Harking, F., Olsen, J.V., Kirchner, A., Lauer, T., Stahlschmidt, M., Talamo, S., Meller, H., Dietl, H., Orschiedt, J., McPherron, S., Krause, J., Meyer, M., Welker, F., Schüler, T., Weiss, M., Hublin, J.-J. 2024. Homo sapiens reached the higher latitudes of Europe by 45,000 years ago. *Nature* 626, 341–346.
- Samantha Brown 2021. Identifying ZooMS Spectra (mammals) using mMass. protocols.io <https://dx.doi.org/10.17504/protocols.io.bzscp6aw>
- Welker, F., Soressi, M., Rendu, W., Hublin, J.-J., & Collins, M. (2015). Using ZooMS to identify fragmentary bone from the Late Middle/Early Upper Palaeolithic sequence of Les Cottés, France. *Journal of Archaeological Science*, 54, 279–286. <https://doi.org/10.1016/j.jas.2014.12.010>

Acknowledgements

We like to thank Anna Wagner and Geoff M. Smith (University of Reading) in their help testing this protocol on archaeological samples.