

Smartphone multiplex microcapillary diagnostics using Cygnus: development and evaluation of rapid serotype-specific NS1 detection with dengue patient samples

Article

Supplemental Material

Creative Commons: Attribution-Share Alike 4.0

Open Access

Needs, S. H. ORCID: <https://orcid.org/0000-0003-3407-9637>,
Sirivisoot, S. ORCID: <https://orcid.org/0000-0002-5809-6936>,
Jegouic, S., Prommool, T., Luangaram, P., Srisawat, C.,
Sriraksa, K., Limpitikul, W., Mairiang, D. ORCID:
<https://orcid.org/0000-0002-2743-0145>, Malasit, P., Avirutnan,
P. ORCID: <https://orcid.org/0000-0003-0053-3045>, Puttikhunt,
C. and Edwards, A. ORCID: <https://orcid.org/0000-0003-2369-989X> (2022) Smartphone multiplex microcapillary diagnostics using Cygnus: development and evaluation of rapid serotype-specific NS1 detection with dengue patient samples. PLOS Neglected Tropical Diseases, 16 (4). e0010266. ISSN 1935-2735 doi: <https://doi.org/10.1371/journal.pntd.0010266>
Available at <https://centaur.reading.ac.uk/104588/>

It is advisable to refer to the publisher's version if you intend to cite from the work. See [Guidance on citing](#).

To link to this article DOI: <http://dx.doi.org/10.1371/journal.pntd.0010266>

Publisher: Public Library of Science

All outputs in CentAUR are protected by Intellectual Property Rights law, including copyright law. Copyright and IPR is retained by the creators or other copyright holders. Terms and conditions for use of this material are defined in the [End User Agreement](#).

www.reading.ac.uk/centaur

CentAUR

Central Archive at the University of Reading

Reading's research outputs online

1 **Supporting Information**

2 **Preparation of MCF multiplex test strips**

3 MCF was produced by Lamina Dielectrics Ltd (Billingshurst, West Sussex, UK). It was
4 made from FEP (Teflon, Dupont, USA) and controlled for size and shape through melt
5 extrusion with air injection to incorporate microcapillaries. MCF was supplied in 100 m reels
6 of flat plastic ribbon containing 10 parallel microcapillaries with an average internal diameter
7 of $206\pm 12.2\ \mu\text{m}$ (Fig 1A). Capture antibodies were diluted in PBS and injected into individual
8 capillaries using a 30G needle (Sigma Aldrich, UK) and incubated at room temperature for 1
9 hour to bind the internal surface by passive adsorption. The antibody solution was removed
10 using a syringe to aspirate air through the capillaries, and replaced with 0.1% (w/v) polyvinyl
11 alcohol (MW 146,000–186,000, > 99% hydrolyzed, Sigma Aldrich, Dorset, UK) in ultra-pure
12 water at room temperature for 2 h (25) to create a hydrophilic internal surface. Coated MCF
13 reels were then blocked with 4% (w/v) BSA (GE Healthcare Life Sciences, Pasching, Austria)
14 overnight at 4 °C. To make the MCF test strip, coated reels were cut into pieces 47 ± 1 mm long
15 and placed into the Cygnus MCF holder cartridges (Fig 1B). The Cygnus holders, strip wells
16 and strip well holder (Figs 1C-E) were 3D printed from PLA (Ooznest, UK) on a Prusa i3 MK3
17 3D printer (Prusa Research, Czech Republic), then washed before each use with detergent,
18 70% ethanol and thorough rinsing in water and air dried; holders were reused multiple times
19 with no loss in performance.

20 Mouse monoclonal anti-DENV NS1 antibody clones which are serotype-specific were
21 chosen from sets previously described and recently modified to be suitable for serotyping-NS1
22 ELISA (9, 17) while the pan-serotype DENV1-4 antibody clone was changed in this study. All
23 clones are listed in Table 1. After monoclonal antibody (mAb) purification, the quality,
24 reactivity, and quantity were examined by gel electrophoresis, western blot and NanoDrop
25 spectrophotometry (Thermo Scientific, IL, USA), respectively. The pan-serotype DENV NS1

26 mAb was biotin labeled (EZ-Link NHS-Biotin; Thermo Scientific) and used as a common
27 detection reagent in all assays. Serotype-specific DENV NS1 mAbs were coated onto the inner
28 surface of MCF at 60 $\mu\text{g/ml}$ for DENV1 and 40 $\mu\text{g/ml}$ for DENV2-4 in duplicate capillaries at
29 room temperature for 1 h. The last two capillaries were used as negative (PBS) and positive
30 (recombinant DENV2 NS1 at 1 $\mu\text{g/ml}$) controls, respectively (Fig 1G).

31 **Enzyme amplified fluorescent ELISA in MCF platform using Cygnus devices**

32 To perform the immunoassays, 60 μl of neat plasma was added into each sample well
33 and a stack of MCF devices were interfaced with the well strip, allowing the samples to flow
34 into each strip and then allowed to incubate for 10 min. The sequence of reagents was similar
35 to conventional microplate ELISA but with significantly shorter incubation steps and is
36 outlined in Fig 1F. Washing buffer (0.05% PBST+0.3 mg/ml of bromophenol blue) in a fresh
37 strip of wells was run through the MCF test strips for 1 min. Incorporation of inert blue dye in
38 the wash buffer allowed visual checking of effective washing and reagent addition.
39 Biotinylated pan-serotype DENV NS1 antibody at 10 $\mu\text{g/ml}$ and a 1:10 dilution of rat serum
40 (to block non-specific background) was added and incubated in the MCF for 5 min followed
41 washing and addition of a dilution 1:1000 of Streptavidin-AP (Sigma Aldrich, Dorset, UK) for
42 5 min. Lastly, a series of three replicate washes was followed by addition of AttoPhos®
43 (Promega, WI, USA) substrate and incubated to allow enzyme amplification of NS1 detected,
44 with images taken after 5, 10 and 15 minutes. The green fluorescence was examined under a
45 blue LED light excitation, imaged through an amber emission filter using blue light
46 transilluminator (IORodeo, USA). The image was recorded by digital camera (Canon S120)
47 and smartphone camera (iPhone 6 and Sony Xperia) within 15 min of substrate addition.
48 Fluorescein at 0, 1, and 5 mM were freshly prepared and used for fluorescent signal
49 normalization by imaging samples alongside reference strips of MCF containing these three

50 dye concentrations, permitting direct comparison of signal intensities between different images
51 and experiments even if camera image properties fluctuated.